



For you to go...



 **KYMCO**
HEALTHCARE



Micro *ForU*

The Micro was designed to be the most transportable, compact and easy to manoeuvre scooter. With puncture proof tyres as standard there is nothing to stop you from going where you want to, without worrying about punctures.



A wireless battery & dismantling system means there are no fiddly plugs or wires to worry about when stripping the scooter down for transporting, in even the smallest of cars. There are even built in handles in the front bumper & rear section to assist handling. With lightweight front & rear sections, loading couldn't be easier. For easy transport, comfort, safety & the ultimate in manoeuvrability, the Micro is *ForU*!



COLOUR OPTIONS

- Blue
- Graphite



MiniLS *ForU*

The new improved Mini LS *ForU* has been designed to offer you exceptional performance while being compact, manoeuvrable & transportable. Incorporating ergonomic design with extra value added features, the Mini LS *ForU* will provide you with a pleasurable and luxurious driving experience - Test drive one today!



This lightweight compact model boasts a whole host of new exciting features. New improved comfort seating & delta handlebars have been added for improved comfort. The latest technology has been utilised for low power LED lighting. There are no plugs or wires to dismantle, so you can strip this scooter in seconds to transport in the boot. Soft-roll puncture proof wheels provide a comfortable & carefree journey. 22ah batteries means more travel range. For range, comfort & portability, the Mini LS is *ForU*!



COLOUR OPTIONS

- Golden Lemon
- Blue
- Graphite
- Flame Orange
- Glossy Black





Super 4

The Super 4 midi scooter is built to handle longer journeys over tougher terrain than other mid size 4 mph mobility scooters. Particular attention has been paid to creating superior styling & comfort. With low profile soft-roll puncture proof tyres there's nothing stopping you! The new high back comfort seat with slide function ensures that you get the most comfortable of journeys!



Super 4 is surprisingly manoeuvrable for a mid size scooter, with a particularly tight turning ability due to a clever steering design developed by KYMCO. With an easy to adjust steering column, comfort seating & extra footboard space, Super 4 provides you with a truly pleasurable & comfortable driving experience. If you need an exceptional all rounder that is comfortable & stylish, Super 4 is ForU



COLOUR OPTIONS

● Blue ● Graphite



Super 8

The Super 8 mid size scooter was designed to offer you the ultimate driving experience for comfort, style & performance. With fully adjustable suspension & a top speed of 8 mph, this scooter offers you the greatest value for money for a compact high performance scooter!



Super 8 is packed with plenty of value added extras to make your driving experience as pleasurable as possible. The fully adjustable suspension system not only allows you to drive safely over rougher ground but can be adjusted for the perfect level of comfort to suit you! With 36 amp batteries and low profile soft-roll puncture proof tyres you won't have to worry about your journey, this scooter will keep you going. For superb performance & great value, the Super 8 is ForU!



COLOUR OPTIONS

● Blue ● Graphite



Midi XLS *ForU*



The New improved Midi XLS now has a fresh new choice of modern colours as well as new features such as adjustable suspension to improve your comfort.

Test drive the new Midi XLS today & feel the difference!



The new Midi XLS *ForU* 8 mph class 3 scooter comes into a category of its own when it comes to comfort & performance. This superb all rounder has many extras included as standard allowing you the comfort & style usually expected on high end Maxi type scooters. The stylish two tone front modesty panel with a convenient built in basket means you can keep your valuables close to hand. Delta hand controls, twin rear view mirrors & a lockable rear basket have been included for your added comfort, safety & security. If you need comfort, style & manoeuvrability on a compact high powered 8 mph scooter, the Midi XLS is the scooter *ForU*!



COLOUR OPTIONS

- Blue ● Graphite
- Flame Orange ● Glossy Black



Maxi XLS *ForU*



When it comes to high performance and absolute comfort, Maxi scooters really do come into a category of their own. Robust enough to tackle the most challenging & toughest environments, you can rest assured that these Maxi models can take you where you want to go, in comfort and style.

The new Maxi XLS is packed with high end features you would expect from this premium luxury scooter. With all round adjustable suspension, these big hitters will handle the most demanding terrain. Delta hand controls, low consumption LED lighting & a rear basket have been included as standard for your convenience. With the Maxi XLS *ForU* you can rest assured that you will be fully equipped to travel in the utmost style, comfort & luxury. If you are looking for a premium scooter, that is fully equipped then the Maxi XLS simply has to be the scooter *ForU*!



COLOUR OPTIONS

- Black ● Blue ● Graphite



Technical Specifications

	Max User Weight	Length	Width	Height	Max Speed	Range	Wheel Size	Heaviest Part (kg)	Total Weight inc Bat (kg)	Battery	Seat Height - deck (mm)	Lights Fitted	Turning radius	Max Slope - Deg
	KG / stone	mm	mm	mm	Kph / mph	Km / Miles	mm (inches)			AH	mm	YES / -	mm	
MicroForU EQ10AA(UK)	114 KG / 18 stone	980 mm	470 mm	830 mm	6.4 Kph / 4 mph	13 Km / 8 Miles	200 mm (7.5")	14.4	42	12AH	330 - 420	-	1110 mm	8
Mini LS ForU EQ20CC(UK)	127 KG / 20 stone	1080 mm	520 mm	860 mm	6.4 Kph / 4 mph	26 Km / 16 Miles	230 mm (9")	19.5	54	22AH	340 - 430	YES	1220 mm	8
Super4 ForU EQ30CB(UK)	127 KG / 20 stone	1220 mm	590 mm	1000 mm	6.4 Kph / 4 mph	32 Km / 20 Miles	255 mm (10")	24	84.4	36AH	420 - 480	YES	1350 mm	10
Super8 ForU EQ35CB(UK)	127 KG / 20 stone	1220 mm	590 mm	1010 mm	12.8 Kph / 8 mph	32 Km / 20 Miles	255 mm (10")	26.1	98.7	36AH	450 - 510	YES	1560 mm	10
Midi XLS ForU EQ35BC(UK)	160 KG / 25 stone	1250 mm	590 mm	1040 mm	12.8 Kph / 8 mph	40 Km / 25 Miles	255 mm (10")	26.1	102	50AH	450 - 510	YES	1560 mm	10
Maxi XLS ForU EQ40BC(UK)	200 KG / 31 stone	1450 mm	670 mm	1090 mm	12.8 Kph / 8 mph	57 Km / 35 Miles	300 mm (12")	52.5	138	80AH	450 - 510	YES	1650 mm	12



Heol Mostyn,
Village Farm Industrial Estate
Pyle, Bridgend CF33 6BJ

t: 01656 670095 | f: 01656 858353

e: sales@kymcohealthcare.co.uk

w: www.kymcohealthcare.co.uk

DEALER STAMP



Some product images include non-standard accessories for illustrative purposes!
Product details are correct at time of print and may change without prior notice!