

MEMPHIS

Rise & Recliner
Dual Motor




Restwell.

*“Memory foam
seat and
dual motor”*

**PETITE
MODEL
NOW
AVAILABLE**

CE

The Memphis Rise Recliner has been designed with comfort in mind. It is the first chair in the Restwell range to feature a memory foam seat. The dual motor mechanism enables the user to find their perfect comfort position. This popular chair is now available in a petite size option.

www.restwellchairs.co.uk

driveTM
www.drivemedical.co.uk

The New Memphis Rise & Recliner is upholstered in our new soft feel Chenille fabric and is available in Berry, Biscuit, Charcoal and Mushroom.



Superb Comfort, Style & Performance...

- > Memory foam seat
- > Soft feel Chenille fabric
- > Stylish T Back design
- > Four-button handset & 2 motors for an independent action
- > Battery back up facility*
- > Ergonomically contoured seat and back
- > Convenient side pocket
- > Dual lift and recline - 3 way mechanism
- > Safety function - key lock only
- > Infinite recline positions
- > Full length chaise
- > Available in standard and petite size options
- > New 2 year warranty**
- > Weight capacity 150kg (23st)

* Batteries not included
 ** Warranty conditions apply, for more information please see user manual

MEMPHIS

CODE: CLR24C (Standard)

CLR34C (Petite)


Max Weight:
150kg (23st)




Height:
112cm (44") - St
104cm (41") - Pe



Width:
84cm (33") - St
77cm (30") - Pe



Depth:
92cm (36") - St
87cm (34") - Pe



Seat Depth:
50cm (19½") - St
46cm (18") - Pe



Seat Height:
48cm (19") - St
46cm (18") - Pe



Seat Width:
53cm (21") - St
48cm (19") - Pe

Colour Options:



Berry



Biscuit



Charcoal



Mushroom

Leaflet Code: LL165

If you would like to view this product please contact your local stockist, details below.



Four button handset